


Air-Lo⁸

Air Alternating Mattress

Bluetooth controlled pressure release
25% better therapy than competing systems

- Powerful yet silent 10 litre motor
- Remote control all functions from a smart phone or tablet
- Quick release heel zone for total pressure off-loading
- Spring loaded hang hooks to secure control unit to bed end



Independent pressure mapping and facility trials demonstrate that an AirLo8 can replace a bi-hourly manual turning regimen*



* May 2021 Independent Pressure mapping Cerebral Palsy Alliance

An optional internal AirLo Hyper-Stretch Coverlet can easily be fitted under all AirLo top covers to ensure maximum patient comfort and insulation

KAREAGE

Call 0470 432 303 or kunal@kareage.com.au
Kareage Pty Ltd www.kareage.com.au

Air-Lo⁸

Specifications



Control Unit:

Dimensions:	W 27 x H 21 x D 10cm
Weight	2.5kg
Mode of Operation	Non-Continuous
Power Rating	240V 50Hz 12VA
Transport Function	Hose End Cap
Air Flow	10 Litres per minute
Auto Startup	✓
Auto Sensor for Patient Weight Profile and Position	✓
Comfort Override Control	✓
Self Diagnostics	✓
Audible/Visual Alarm	✓
Mute Function	✓
Static Mode	Timeout - 60 minutes
Care Mode	Timeout - 20 minutes
Incline mode with Auto Tilt Sensor	✓
Auto Anti Tamper - Control Panel Unlock Function	✓
CE Certification and C-Tick Listed N27882	✓

Warranty

3 Years Pump
3 Year Soft Goods

ARTG: 175810

Mattress Replacement

Dimensions Standard	L 200 x W 90 x H 23cm
Dimensions King	L 200 x W 107 x H 23cm
Alternating Cell Type	Safety Cell in Cell
Overall Cells Count	20 - (3 static Head Cells included)
Quick release valves for heel zone cells	6
Cell Material	Pure TPU
Top Cover	Dahlia Bi-Elastic 2 way stretch - welded
Base/s	Nylon TPU welded
Cell Cycle	1 in 2
Alternating Mode Default	10 minutes per alternation
Infection Control	Welded Top Cover
Positioning Handles	2 aside
Base Attachment	Side Straps

Optional Hyperstretch Coverlet

Dimensions	L 200 x W 90/107 x H 1cm
Attachment	Base Surround Zip
Cover Material	High Performance Micromesh
Foam	Dunlop Marathon Hypersoft HS18-35 with Ultrafresh